MICROBIX

REDx™ FLOQ®

REDx™ FLOQ® SARS-CoV-2/Flu A&B/RSV Swab Positive Control

Microbix Biosystems Inc.
265 Watline Avenue
Mississauga, Ontario,
Canada, L4Z 1P3

Cat#: RED-S-19-M1

FLOQ® is a trademark of Copan Italia Spa

OBELIS S.A.
Bd. General Wahis, 53
1030 Brussels
Belgium
+ 32 2 732 59 54
mail@obelis.com

About this package insert

Thank you for your interest in this REDx™ quality control product.

This package insert consists of two pages.

- The first page contains the product name and an explanation of the symbols used on the labeling.
- The second page contains the complete package insert text.

If the package insert you view or print does not contain two pages, or if you experience any problems, please email us at customer.service@microbix.com.

By phone: US customers call +1-800-794-6694; International customers call collect +1-905-361-8910.

A printed package insert will be sent to you upon request.

P/N RED-S-19-M1.5R0

Explanation of symbols used in Microbix product labeling

Upper limit
of temperature

Temperature
limitation

In Vitro Diagnostic
Medical Device

Single-use only

Positive control

Use By

"Caution, consult
accompanying documents"

Catalogue
number

Batch code

Manufacturer

WARNING: THESE REAGENTS MUST NOT BE SUBSTITUTED FOR THE MANDATORY POSITIVE AND NEGATIVE SAMPLE REAGENTS PROVIDED WITH MANUFACTURED TEST KITS.

REDx™ FLOQ® SARS-CoV-2/Flu A&B/RSV Swab Positive Control

FOR IVD USE.

INTENDED USE

REDx™ FLOQ® SARS-CoV-2/Flu A&B/RSV Swab Positive Control is a desiccated, unassayed control intended to monitor laboratory testing performance, procedures, and workflow with nucleic acid tests (NATs) that detect SARS-CoV-2, Influenza A, Influenza B, and/or Respiratory Syncytial Virus in human nasopharyngeal, oropharyngeal, nasal mid-turbinate, anterior nares and lower respiratory samples, collected on swabs.

PRODUCT DESCRIPTION

REDx™ FLOQ® SARS-CoV-2/Flu A&B/RSV Swab Positive Control is formulated with inactivated whole-genome SARS-CoV-2 cDNA and inactivated native Influenza A, Influenza B, and Respiratory Syncytial Virus. REDx™ FLOQ® SARS-CoV-2/Flu A&B/RSV Swab Positive Control can be utilized as an external sample to monitor the process of respiratory virus nucleic acid detection assays, including sample extraction and purification, amplification, and detection¹.

REDx™ FLOQ® SARS-CoV-2/Flu A&B/RSV Swab Positive Control does not have an assigned value ("unassayed"). Laboratories are required to establish an acceptance range for each lot of REDx™ FLOQ® SARS-CoV-2/Flu A&B/RSV Swab Positive Control with all assay procedures that the control is intended to be used with, prior to routine use in the laboratory^{2,3}.

PRINCIPLES OF THE PROCEDURE

REDx™ FLOQ® SARS-CoV-2/Flu A&B/RSV Swab Positive Control is designed as an external independent sample for use with laboratory testing of SARS-CoV-2, Influenza A, Influenza B, and/or Respiratory Syncytial Virus nucleic acid targets, according to ISO 15189 and CLIA regulations.

REAGENTS

Cat. No RED-S-19-M1; 1 swab containing inactivated Influenza A, Influenza B, Respiratory Syncytial Virus, and whole-genome SARS-CoV-2 cDNA.

LIMITATIONS OF THE PROCEDURE

REDx™ CONTROLS MUST NOT BE SUBSTITUTED FOR THE POSITIVE AND NEGATIVE CONTROL REAGENTS PROVIDED WITH THE MANUFACTURED TEST KITS.

TEST PROCEDURES and INTERPRETATION OF RESULTS provided by manufacturers of test kits must be followed closely.

Deviations from procedures recommended by test kit manufacturers may produce unreliable results. REDx™ FLOQ® SARS-CoV-2/Flu A&B/RSV Swab Positive Control DOES NOT HAVE AN ASSIGNED VALUE and may not be suitable for use with all respiratory viral test kits and procedures. Procedures for implementing a quality assurance program and monitoring test performance on a routine basis must be established by each individual laboratory. Each laboratory should establish its own range of acceptable values^{2,3}.

Samples are not calibrators and should not be used for assay calibration.

REDx™ FLOQ® SARS-CoV-2/Flu A&B/RSV Swab Positive Control is recommended for use with NATs only.

Adverse shipping and storage conditions or use of outdated samples may produce erroneous results. REDx™ FLOQ® SARS-CoV-2/Flu A&B/RSV Swab Positive Control evaluates only the PCR steps of the RT-PCR methods for detection of SARS-CoV-2.

REDx™ FLOQ® SARS-CoV-2/Flu A&B/RSV Swab Positive Control might not be suitable for NATs without an extraction step.

WARNINGS AND PRECAUTIONS

For IVD use.

For Professional and Trained Laboratory Personnel Use Only

Safety Precautions

1. Raw material used for REDx™ FLOQ® SARS-CoV-2/Flu A&B/RSV Swab Positive Control preparation is inactivated.
2. Use Centers for Disease Control and Prevention (CDC) recommended universal precautions for handling the samples and human specimens⁴.
3. REDx™ FLOQ® SARS-CoV-2/Flu A&B/RSV Swab Positive Control must be disposed of by following RCRA ID#D001 guidelines for ignitable waste⁵.
4. Keep REDx™ FLOQ® SARS-CoV-2/Flu A&B/RSV Swab Positive Control pouch closed when not in use.

HANDLING PRECAUTIONS

1. Do not use controls beyond the expiration date.
2. Avoid contamination of controls when opening the swab pouches.

STORAGE INSTRUCTIONS

Store REDx™ FLOQ® SARS-CoV-2/Flu A&B/RSV Swab Positive Control at 2-30°C until use. Once opened, REDx™ FLOQ® SARS-CoV-2/Flu A&B/RSV Swab Positive Control should not be reused.

MATERIALS PROVIDED

REDx™ FLOQ® SARS-CoV-2/Flu A&B/RSV Swab Positive Control – 1 swab

MATERIALS REQUIRED, BUT NOT PROVIDED

Refer to the instructions supplied by the manufacturer of the test kit to be used.

PROCEDURE

When including the REDx™ FLOQ® SARS-CoV-2/Flu A&B/RSV Swab Positive Control in a test run, the exact same procedure for unknown specimens collected on a swab must be used. Refer to the manufacturer's supplied instructions for use provided with the respiratory viral test kit.

1. Elute the REDx™ FLOQ® SARS-CoV-2/Flu A&B/RSV Swab Positive Control by referring to the preferred technique and volumes described in the assay instructions for use (usually 1-3 mL).
2. Use 100-1000 µL from the eluted swab REDx™ FLOQ® SARS-CoV-2/Flu A&B/RSV Swab Positive Control for the nucleic acid extraction step.
3. After extraction, proceed with the NAT by using the eluted nucleic acid test volume specified in the assay instruction for use (usually 5-20 µL from the eluted purified nucleic acid volume).

NOTE: Controls must NOT be substituted for the internal kit positive and negative controls.

As the REDx™ FLOQ® SARS-CoV-2/Flu A&B/RSV Swab Positive Control does not have an assigned value, the laboratory must establish an acceptance range for each lot of the REDx™ FLOQ® SARS-CoV-2/Flu A&B/RSV Swab Positive Control.

TROUBLESHOOTING

When REDx™ FLOQ® SARS-CoV-2/Flu A&B/RSV Swab Positive Control results are outside of the established laboratory acceptance range for internal controls, it may be an indication of unsatisfactory test performance.

Possible sources of error include deterioration of test kit reagents, operator error, faulty performance of equipment, or contamination of reagents; internal laboratory procedures should be followed.

REFERENCES

1. *Accurate Results in the Clinical Laboratory 2013*, ISBN: 978-0-12-415783-5
2. Kinns H, Pitkin S, Housley D, et al. *J Clin Pathol* 2013;66:1027-1032.
3. *Statistical Quality Sample for Quantitative Measurements: Principles and Definitions; Approved Guideline— Second Edition*. NCCLS document C24-A2, 1999.
4. *CDC Recommendations for prevention of HIV transmission in health care settings*. *MMWR* 36 (suppl. 2), 1987.
5. *Treatment standards for hazardous waste; 40 CFR 268.40 Subpart D*. D001: Ignitable characteristics of waste.

For assistance, contact Microbix Technical Support at +1-905-361-8910.

P/N RED-S-19-M1.5R0
17 November, 2021

CERTIFICATE OF ANALYSIS

REDx™ FLOQ® SARS-CoV-2/FluA&B/RSV Swab Positive Control

Product Information

Product Name	REDx™ FLOQ® SARS-CoV-2/FluA&B/RSV Swab Positive Control
Catalogue Number	RED-S-19-M1
Strain/Cell line	<i>Influenza A/H1N1 Subtype; Influenza B/Hong Kong; Human respiratory syncytial virus, long strain; SARS-CoV-2 whole genome cDNA (NC_045512)</i>
Lot number	19M10013A8
Analyte Value	Unassayed. Value not assigned by the manufacturer
Product Description	Unassayed qualitative positive control for the detection of SARS-CoV-2, Influenza A, Influenza B and RSV
Elution Volume	1-3 mL
Storage conditions	2-30°C
Expiry Date	08 Jun 2027
Intended Use	For <i>In Vitro</i> Diagnostic Use

Quality Control Information

Test method	Xpert® Xpress SARS-CoV-2/Flu/RSV
Test Result	Observed Result
SARS-CoV-2	Positive
Influenza A	Positive
Influenza B	Positive
RSV	Positive
Inactivation	Inactivated

This product is an unassayed control. The reported QC information is not intended to represent product specifications in any commercial or lab developed test. The laboratory should establish its own analyte value and ranges.

Microbix Biosystems Inc.
 265 Watline Avenue
 Mississauga, Ontario, L4Z 1P3
 Obelis S.A.
 Bd. General Wahis, 53 1030 Brussels
 Belgium

Information

Please ensure that your institution is consulted before recycling an item.

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RÉFÉRENCES

1. *Accurate Results in the Clinical Laboratory 2013*, ISBN: 978-0-12-12-415783-5
2. *Kinns H, Pitkin S, Housley D, et al. J. Clin Pathol 2013;66:1027-1032.*
3. *Statistical Quality Sample for Quantitative Measurements: Principles and Definitions; Approved Guideline— Second Edition. NCCLS document C24-A2, 1999.*
4. *CDC Recommendations for prevention of HIV transmission in health care settings. MMWR 36 (suppl. 2), 1987.*
5. *Treatment standards for hazardous waste; 40 CFR 268.40 Subpart D. D001: Ignitable characteristics of waste.*

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REFERENZEN

- Accurate Results in the Clinical Laboratory 2013*, ISBN: 978-0-12-415783-5
- Kinns H, Pitkin S, Housley D, et al. *J Clin Pathol* 2013;66:1027–1032.
- Statistical Quality Sample for Quantitative Measurements: Principles and Definitions; Approved Guideline— Second Edition*. NCCLS document C24-A2, 1999.
- CDC Recommendations for prevention of HIV transmission in health care settings*. *MMWR* 36 (suppl. 2), 1987.
- Treatment standards for hazardous waste; 40 CFR 268.40 Subpart D. D001: Ignitable characteristics of waste*.

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REFERENCIAS

1. *Accurate Results in the Clinical Laboratory 2013*, ISBN: 978-0-12-415783-8
2. Kinns H, Pitkin S, Housley D, et al. *J Clin Pathol* 2013;66:1027-1032.
3. *Statistical Quality Sample for Quantitative Measurements: Principles and Definitions; Approved Guideline- Second Edition*. NCCLS document C24-A2, 1999.
4. *CDC Recommendations for prevention of HIV transmission in health care settings*. *MMWR* 36 (suppl. 2), 1987.
5. *Treatment standards for hazardous waste; 40 CFR 268.40 Subpart D. D001: Ignitable characteristics of waste*.

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- Accurate Results in the Clinical Laboratory 2013*, ISBN: 978-0-12-415783-5
- Kinns H, Pitkin S, Housley D, et al. J Clin Pathol 2013;66:1027-1032.*
- Statistical Quality Sample for Quantitative Measurements: Principles and Definitions; Approved Guideline- Second Edition. NCCLS document C24-A2, 1999.*
- CDC Recommendations for prevention of HIV transmission in health care settings. MMWR 36 (suppl. 2), 1987.*
- Treatment standards for hazardous waste; 40 CFR 268.40 Subpart D. D001: Ignitable characteristics of waste.*

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

Information  Please ensure that your institution is consulted before recycling an item.

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Quality Assurance Signature: *Jenn Lee*

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BRONNEN

1. *Accurate Results in the Clinical Laboratory 2013*, ISBN: 978-0-12-415783-5
2. *Kinns H, Pitkin S, Housley D, et al. J Clin Pathol 2013;66:1027-1032.*
3. *Statistical Quality Control for Quantitative Measurements: Principles and Definitions; Approved Guideline— Second Edition. NCCLS document C24-A2, 1999.*
4. *CDC Recommendations for prevention of HIV transmission in health care settings. MMWR 36 (suppl. 2), 1987.*
5. *Treatment standards for hazardous waste; 40 CFR 268.40 Subpart D. D001: Ignitable characteristics of waste.*

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R-Biopharm AG