MICROBIX

REDxTM HPV 18 Positive Control

Microbix Biosystems Inc.
265 Watline Avenue,
Mississauga Ontario,
Canada L4Z 1P3

Cat#: RED-62-18

About this package insert

Thank you for your interest in this REDxTM Controls product.

This package insert consists of two pages.

- The first page contains the product name and an explanation of the symbols used on the labeling.
- The second page contains the complete package insert text.

If the package insert you view or print does not contain two pages, or if you experience any problems, please email us at customer.service@microbix.com. By phone: US customers call +1-800-794-6694; International customers call collect +1-905-361-8910.

A printed package insert will be sent to you upon request.

OBELIS S.A
Bd. General Wahis, 53
1030 Brussels
Belgium
+ 32 2 732 59 54
mail@obelis.com

P/N RED-62-18.5R2

Explanation of symbols used in Microbix product labeling

Upper limit
of temperature

Temperature
limitation

Highly flammable

In Vitro Diagnostic
Medical Device

Biological risks

Use By

"Caution, consult
accompanying documents"

Toxic by inhalation, in contact
with skin and if swallowed

Positive control

Catalogue
number

Authorized Representative in
the European Community

Single-use only

Batch code

Manufacturer

WARNING: THESE REAGENTS MUST NOT BE SUBSTITUTED FOR THE MANDATORY POSITIVE AND NEGATIVE SAMPLE REAGENTS PROVIDED WITH MANUFACTURED TEST KITS.

REDx™ HPV 18 Positive Control

FOR IVD USE.

INTENDED USE

REDx™ HPV 18 Positive Control is an unassayed control intended to monitor laboratory testing performance, procedures, and workflow with nucleic acid tests (NATs) that detect high-risk Human Papillomavirus (HPV) type 18 DNA and/or RNA in human cervical, anal, penile, or oropharyngeal samples collected in an alcohol-based transport medium.

PRODUCT DESCRIPTION

REDx™ HPV 18 Positive Control is formulated with whole human epithelial cells with integrated and episomal HPV 18 DNA and/or RNA presence. REDx™ HPV 18 Positive Control can be utilized as an external sample to monitor the process of HPV nucleic acid detection assays, including sample extraction and purification, amplification, and detection¹.

REDx™ HPV 18 Positive Control does not have an assigned value ("unassayed"). Laboratories are required to establish an acceptance range for each lot of REDx™ HPV 18 Positive Control with all assay procedures that the control is intended to be used with, prior to routine use in the laboratory^{2,3}.

PRINCIPLES OF THE PROCEDURE

REDx™ HPV 18 Positive Control is designed as an external independent sample for use with laboratory testing of HPV 18 DNA and/or RNA targets, according to ISO 15189 and CLIA regulations.

REAGENTS

Cat. No RED-62-18; 1 vial formulated with cells containing human epithelial cells with integrated and episomal HPV 18 DNA and/or HPV 18 RNA presence.

NOTE: REDx™ HPV 18 Positive Control contains a suspension of cells in an alcohol solution and may therefore, exhibit slight cloudiness.

LIMITATIONS OF THE PROCEDURE

REDx™ CONTROLS MUST NOT BE SUBSTITUTED FOR THE POSITIVE AND NEGATIVE CONTROL REAGENTS PROVIDED WITH THE MANUFACTURED TEST KITS.

TEST PROCEDURES and INTERPRETATION OF RESULTS provided by manufacturers of test kits must be followed closely.

Deviations from procedures recommended by test kit manufacturers may produce unreliable results.

REDx™ HPV 18 Positive Control DOES NOT HAVE AN ASSIGNED VALUE and may not be suitable for use with all HPV test kits and procedures. Specific levels of reactivity will vary among different manufacturers' assays, different procedures, and different laboratories. Procedures for implementing a quality assurance program and monitoring test performance on a routine basis must be established by each individual laboratory. Each laboratory should establish its own range of acceptable values^{2,3}.

Samples are not calibrators and should not be used for assay calibration.

REDx™ HPV 18 Positive Control is recommended for use with NATs only.

Adverse shipping and storage conditions or use of outdated samples may produce erroneous results.

REDx™ HPV 18 Positive Control might not be suitable for NATs without an extraction step.

WARNINGS AND PRECAUTIONS

For IVD use.

For Professional and Trained Laboratory Personnel Use Only

WARNING: Handle controls as if they are capable of transmitting infectious agents. REDx™ HPV 18 Positive Control is manufactured from HPV infected cells and other non-infected cells that are grown in tissue culture and preserved in a buffered solution.

Safety Precautions

1. Raw material used for REDx™ HPV 18 Positive Control preparation is inactivated.
2. Use Centers for Disease Control and Prevention (CDC) recommended universal precautions for handling the samples and human specimens⁴.
3. REDx™ HPV 18 Positive Control must be disposed of by following RCRA ID#D001 guidelines for ignitable waste⁵.
4. Do not pipette by mouth; do not smoke, eat or drink in areas where specimens are being handled. Clean any spillage by immediately wiping up with 0.5% sodium hypochlorite solution. Dispose of all specimens, controls, and materials used in testing as though they contain infectious agents.
5. Keep REDx™ HPV 18 Positive Control closed when not in use. Avoid direct inhalation of the solution and use with ventilation.

HANDLING PRECAUTIONS

1. Do not use controls beyond the expiration date.
2. Avoid contamination of controls when opening and closing the vials.
3. REDx™ HPV 18 Positive Control contains a FLAMMABLE liquid; keep away from all sources of ignition.

STORAGE INSTRUCTIONS

Store REDx™ HPV 18 Positive Control at 2-8°C until use.

Once opened, REDx™ HPV 18 Positive Control should not be reused. Store vials upright to prevent leakage.

MATERIALS PROVIDED

REDx™ HPV 18 Positive Control – 1mL vial x1

MATERIALS REQUIRED, BUT NOT PROVIDED

Refer to the instructions supplied by the manufacturer of the test kit to be used.

PROCEDURE

When including the REDx™ HPV 18 Positive Control in a test run, the exact same procedure for unknown specimens collected in an alcohol-based transport medium must be used.

Before use, thoroughly mix the REDx™ HPV 18 Positive Control by aspirating 0.5 mL with a P1000 micropipette and pulse vortexing for at least 30 seconds.

Refer to the manufacturer's supplied instructions for use provided with the HPV test kit.

NOTE: Controls must NOT be substituted for the internal kit positive and negative controls. REDx™ HPV 18 Positive Control is an independent control intended to monitor laboratory testing performance. However, in case of any serious incident that has occurred in relation to the product, report the event to Microbix and the competent authority of the country in which the event occurred.

Levels of reactivity of REDx™ HPV 18 Positive Control may vary with different manufacturer's test kits and different kit lots. As REDx™ HPV 18 Positive Control does not have an assigned value, the laboratory must establish an acceptance range for each lot of REDx™ HPV 18 Positive Control.

TROUBLESHOOTING

When REDx™ HPV 18 Positive Control results are outside of the established laboratory acceptance range for internal controls, it may be an indication of unsatisfactory test performance. Possible sources of error include deterioration of test kit reagents, operator error, faulty performance of equipment, or contamination of reagents; internal laboratory procedures should be followed.

REFERENCES

1. *Accurate Results in the Clinical Laboratory 2013*, ISBN: 978-0-12-415783-5
2. Kinns H, Pitkin S, Housley D, et al. *J Clin Pathol* 2013;66:1027–1032.
3. *Statistical Quality Sample for Quantitative Measurements: Principles and Definitions; Approved Guideline— Second Edition. NCCLS document C24-A2, 1999.*
4. *CDC Recommendations for prevention of HIV transmission in health care settings. MMWR* 36 (suppl. 2), 1987.
5. *Treatment standards for hazardous waste; 40 CFR 268.40 Subpart D. D001: Ignitable characteristics of waste.*

For assistance, contact Microbix Technical Support at +1-905-361-8910.

P/N RED-62-18.5R2
08 October, 2024

265 Watline Ave.
Mississauga, ON
Canada L4Z 1P3
Tel. 905-361-8910
Fax 905-361-8911
microbix.com

CERTIFICATE OF ANALYSIS

REDx™ HPV 18 Positive Control

Product Information

Product Name	REDx™ HPV 18 Positive Control
Catalogue Number	RED-62-18
Strain	HPV type 18
Lot number	62180013A1
Analyte Value	Unassayed. Value not assigned by the manufacturer
Product Description	Unassayed qualitative positive control for the detection of HPV 18. Contains human DNA
Fill Volume	1 mL
Buffer	Aqueous solution containing methanol as preservative
Storage conditions	2-8°C
Expiry Date	30 Apr 2028
Intended Use	For In Vitro Diagnostic Use.

Quality Control Information


Product name	REDx™ HPV 18 Positive Control
Test method	Cepheid GeneXpert HPV test
Test Result	Observed Result
HPV 18/45	Positive
hDNA	Positive
HPV 16	Negative
Inactivation	Inactivated

This product is an unassayed control. The reported QC information is not intended to represent product specifications in any commercial or lab developed test. The laboratory should establish its own analyte value and ranges.

Information Please ensure that your institution is consulted before recycling an item.

Microbix Biosystems Inc.
265 Watline Avenue
Mississauga, Ontario, L4Z 1P3
Obelis S.A.
Bd. General Wahis, 53 1030 Brussels
Belgium

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RÉFÉRENCES

1. *Accurate Results in the Clinical Laboratory 2013*, ISBN: 978-0-12-12-415783-5
2. *Kinns H, Pitkin S, Housley D, et al. J. Clin Pathol 2013;66:1027-1032.*
3. *Statistical Quality Sample for Quantitative Measurements: Principles and Definitions; Approved Guideline— Second Edition. NCCLS document C24-A2, 1999.*
4. *CDC Recommendations for prevention of HIV transmission in health care settings. MMWR 36 (suppl. 2), 1987.*
5. *Treatment standards for hazardous waste; 40 CFR 268.40 Subpart D. D001: Ignitable characteristics of waste.*

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REFERENCIAS

1. *Accurate Results in the Clinical Laboratory 2013*, ISBN: 978-0-12-415783-5
2. *CDC Recommendations for prevention of HIV transmission in health care settings*. *MMWR* 36 (suppl. 2), 1987.
3. *Treatment standards for hazardous waste; 40 CFR 268.40 Subpart D. D001: Ignitable characteristics of waste*
4. Kinns H, Pitkin S, Housley D, et al. *J Clin Pathol* 2013;66:1027-1032.
5. *Statistical Quality Sample for Quantitative Measurements: Principles and Definitions; Approved Guideline— Second Edition*. NCCLS document C24-A2, 1999.

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
Microbix Biosystems Inc.
265 Watline Avenue
Mississauga, Ontario, L4Z 1P3
Obelis S.A.
Bd. General Wahis, 53 1030 Brussels
Belgium
Please ensure that your institution is consulted before recycling an item.

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
REFERENZEN

- Accurate Results in the Clinical Laboratory 2013*, ISBN: 978-0-12-415783-5
- CDC Recommendations for prevention of HIV transmission in health care settings*. *MMWR* 36 (suppl. 2), 1987.
- Treatment standards for hazardous waste; 40 CFR 268.40 Subpart D. D001: Ignitable characteristics of waste*.
- Kinns H, Pitkin S, Housley D, et al. J Clin Pathol 2013;66:1027-1032*.
- Statistical Quality Sample for Quantitative Measurements: Principles and Definitions; Approved Guideline— Second Edition. NCCLS document C24-A2, 1999*.

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1. *Accurate Results in the Clinical Laboratory 2013*, ISBN: 978-0-12-415783-5
2. *CDC Recommendations for prevention of HIV transmission in health care settings*. MMWR 36 (suppl. 2), 1987
3. *Treatment standards for hazardous waste; 40 CFR 268.40 Subpart D. D001: Ignitable characteristics of waste.*
4. *Kinns H, Pitkin S, Housley D, et al. J Clin Pathol 2013;66:1027-1032.*
5. *Statistical Quality Sample for Quantitative Measurements: Principles and Definitions; Approved Guideline— Second Edition. NCCLS document C24-A2, 1999.*

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BRONNEN

1. *Accurate Results in the Clinical Laboratory 2013*, ISBN: 978-0-12-415783-5
2. *CDC Recommendations for prevention of HIV transmission in health care settings*. MMWR 36 (suppl. 2), 1987.
3. *Treatment standards for hazardous waste; 40 CFR 268.40 Subpart D. D001: Ignitable characteristics of waste*.
4. Kinns H, Pitkin S, Housley D, et al. *J Clin Pathol* 2013;66:1027–1032.
5. *Statistical Quality Control for Quantitative Measurements: Principles and Definitions; Approved Guideline— Second Edition*. NCCLS document C24-A2, 1999.

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R-Biopharm AG