REDx™FLOQ® Respiratory Swab Negative Control

Microbix Biosystems Inc.
265 Watline Avenue
Mississauga, Ontario, Canada, L4Z 1P3

Cat#: RED-S-99-01

FLOQ® is a trademark of Copan Italia Spa

OBELIS S.A
Bd. General Wahis, 53
1030 Brussels
Belgium
+ 32 2 732 59 54
mail@obelis.com

About this package insert

Thank you for your interest in this REDx™ quality control product. This package insert consists of two pages.

- The first page contains the product name and an explanation of the symbols used on the labeling.
- The second page contains the complete package insert text.

If the package insert you view or print does not contain two pages, or if you experience any problems, please email us at customer.service@microbix.com.

By phone: US customers call +1-800-794-6694; International customers call collect

+1-905-361-8910.

A printed package insert will be sent to you upon request.

P/N RED-S-99-01.5R1

Explanation of symbols used in Microbix product labeling

Upper limit
of temperature

Temperature
limitation

Authorized Representative in
the European Community

In Vitro Diagnostic
Medical Device

Negative control

Use By

"Caution, consult
accompanying documents"

Single-use only

Catalogue
number

Batch code

Manufacturer

WARNING: THESE REAGENTS MUST NOT BE SUBSTITUTED FOR THE MANDATORY POSITIVE AND NEGATIVE SAMPLE REAGENTS PROVIDED WITH MANUFACTURED TEST KITS.

REDx™FLOQ® Respiratory Swab Negative Control

FOR IVD USE.

INTENDED USE

REDx™FLOQ® Respiratory Swab Negative Control is a desiccated, unassayed control intended to monitor laboratory testing performance, procedures, and workflow with diagnostic tests that detect respiratory pathogens in human nasopharyngeal, oropharyngeal, nasal mid-turbinate and anterior nares samples, collected on swabs.

PRODUCT DESCRIPTION

REDx™FLOQ® Respiratory Swab Negative Control can be utilized as an external sample to monitor the process of respiratory pathogen diagnostic assays.¹

REDx™FLOQ® Respiratory Swab Negative Control does not have an assigned value (“unassayed”). It is required that each laboratory establish an acceptance range for each lot of Respiratory Swab Negative Control with each assay procedure, with which it is intended to be run, prior to routine use in the laboratory.^{3,4}

PRINCIPLES OF THE PROCEDURE

REDx™FLOQ® Respiratory Swab Negative Control is designed as an external independent sample for use with laboratory testing of respiratory pathogens according to ISO15189 and CLIA regulations. REDx™FLOQ® Respiratory Swab Negative Control is manufactured with a proprietary matrix and excipients.

REAGENTS

Cat. No RED-S-99-01 1 swab containing proprietary matrix and excipients.

LIMITATIONS OF THE PROCEDURE

REDx™ SAMPLES MUST NOT BE SUBSTITUTED FOR THE POSITIVE AND NEGATIVE SAMPLE REAGENTS PROVIDED WITH MANUFACTURED TEST KITS.

TEST PROCEDURES and INTERPRETATION OF RESULTS provided by manufacturers of test kits must be followed closely.

Deviations from procedures recommended by test kit manufacturers may produce unreliable results.

REDx™ FLOQ® Respiratory Swab Negative Control DOES NOT HAVE AN ASSIGNED VALUE and may not be suitable for use with all respiratory pathogens test kits and procedures. Specific levels of reactivity will vary among different manufacturers’ assays, different procedures and different laboratories. Procedures for implementing a quality assurance program and monitoring test performance on a routine basis must be established by each individual laboratory. Each laboratory is required to establish its own range of acceptable values.^{3,4}

Samples are not calibrators and should not be used for assay calibration.

Adverse shipping and storage conditions or use of outdated Samples may produce erroneous results.

WARNINGS AND PRECAUTIONS

For IVD use.

For Professional and Trained Laboratory Personnel Use Only

Safety Precautions

1. Use Centers for Disease Control and Prevention (CDC) recommended universal precautions for handling the samples and human specimens.²
2. Keep REDx™FLOQ® Respiratory Swab Negative Control pouch sealed prior to use.

Handling Precautions

1. Do not use samples beyond the expiration date.
2. Avoid contamination of samples when opening the swab pouches.

STORAGE INSTRUCTIONS

Store REDx™FLOQ® Respiratory Swab Negative Control at 2-30°C until use.

Once opened REDx™FLOQ® SARS-CoV-2 Swab Negative Control should not be reused.

MATERIALS PROVIDED

REDx™FLOQ® Respiratory Swab Negative Control – 1 swab

MATERIALS REQUIRED, BUT NOT PROVIDED

Refer to the instructions supplied by manufacturer of the test kit to be used.

procedure for unknown specimens collected in on a swab must be used. Refer to the manufacturer supplied instructions for use provided with the test kit.

1. Elute the REDx™FLOQ® Respiratory Swab Negative Control by referring to the preferred technique and volumes described in the assay instructions for use (usually 0.2-0.7 mL).
2. After elution, proceed with detection step by using the test volume specified in the assay instruction for use (usually 50-200 µL from the eluted swab volume).

NOTE: Samples must NOT be substituted for the internal kit positive and negative Samples. REDx™FLOQ® Respiratory Swab Negative Control is an independent control intended to monitor laboratory testing performance. However, in case of any serious incident that has occurred in relation to the product, report the event to Microbix and the competent authority of the country in which the event occurred.

As REDx™FLOQ® Respiratory Swab Negative Control does not have an assigned value, the laboratory is required to establish an acceptance range for each lot of REDx™FLOQ® Respiratory Swab Negative Control.

TROUBLESHOOTING

When results REDx™FLOQ® Respiratory Swab Negative Control are outside of the established laboratory acceptance range, it may be an indication of unsatisfactory test performance.

Possible sources of error include: deterioration of test kit reagents, operator error, faulty performance of equipment, or contamination of reagents; internal laboratory procedures should be followed.

REFERENCES

1. *Accurate Results in the Clinical Laboratory 2013*, ISBN: 978-0-12-415783-5
2. *CDC Recommendations for prevention of HIV transmission in health care settings. MMWR 36 (suppl. 2), 1987.*
3. *Kinns H, Pitkin S, Housley D, et al. J Clin Pathol 2013;66:1027–1032.*
4. *Statistical Quality Sample for Quantitative Measurements: Principles and Definitions; Approved Guideline– Second Edition. NCCLS document C24-A2, 1999.*

For assistance, contact Microbix Technical Support at +1-905-361-8910.

PROCEDURE

When including the REDx™FLOQ® Respiratory Swab Negative Control in a test run, the exact same

P/N RED-S-99-01.5R1

08 October 2024

265 Watline Ave.
Mississauga, ON
Canada L4Z 1P3
Tel. 905-361-8910
Fax 905-361-8911
microbix.com

CERTIFICATE OF ANALYSIS

REDx™FLOQ® Respiratory Swab Negative Control

Product Information

Product Name	REDx™FLOQ® Respiratory Swab Negative Control
Catalogue Number	RED-S-99-01
Lot number	99010038A17
Analyte Value	Unassayed. Value not assigned by the manufacturer
Product Description	Unassayed qualitative negative control for the detection of SARS-CoV-2.
Elution Volume	0.2-0.5 mL
Storage conditions	2-30°C
Expiry Date	15 Sep 2027
Intended Use	For In Vitro Diagnostic Use.

Quality Control Information

Test method	Sofia SARS Antigen FIA
Test Result	Observed Result
SARS-CoV-2	Negative

This product is an unassayed control. The reported QC information is not intended to represent product specifications in any commercial or lab developed test. The laboratory should establish its own analyte value and ranges.

Microbix Biosystems Inc.
265 Watline Avenue
Mississauga, Ontario, L4Z 1P3
Obelis S.A.
Bd. General Wahis, 53 1030|Brussels
Belgium

Information

Please ensure that your institution is consulted before recycling an item.

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RÉFÉRENCES

1. *Accurate Results in the Clinical Laboratory 2013*, ISBN: 978-0-12-12-415783-5
2. *CDC Recommendations for prevention of HIV transmission in health care settings. MMWR 36 (suppl. 2), 1987.*
3. *Kinns H, Pitkin S, Housley D, et al. J Clin Pathol 2013;66:1027–1032.*
4. *Statistical Quality Sample for Quantitative Measurements: Principles and Definitions; Approved Guideline—Second Edition. NCCLS document C24-A2, 1999.*

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
Quality Assurance Signature: *jenn lee*

Date: 22 Apr 2026

MICROBIX

Microbix Biosystems Inc.
265 Watline Avenue
Mississauga, Ontario, Canada, L4Z 1P3

REDTMFLOQ[®]

REDxTMFLOQ[®] Respiratory Swab Negative Control

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R-Biopharm AG