MICROBIX

REDx™ STI Negative Control

Microbix Biosystems Inc.
265 Watline Avenue
Mississauga, Ontario, Canada, L4Z 1P3

Cat#: RED-99-M1

OBELIS S.A
Bd. General Wahis, 53
1030 Brussels
Belgium
+ 32 2 732 59 54
mail@obelis.com

About this package insert

Thank you for your interest in this REDx™ quality control product.
This package insert consists of two pages.

- The first page contains the product name and an explanation of the symbols used on the labeling.
- The second page contains the complete package insert text.

If the package insert you view or print does not contain two pages, or if you experience any problems, please email us at customer.service@microbix.com.
By phone: US customers call +1-800-794-6694; International customers call collect +1-905-361-8910.

A printed package insert will be sent to you upon request.

P/N RED-99-M1.5R1

Explanation of symbols used in Microbix product labeling

Upper limit of temperature

Temperature limitation

Highly flammable

In Vitro Diagnostic Medical Device

Negative control

Use By

"Caution, consult accompanying documents"

Toxic by inhalation, in contact with skin and if swallowed

Manufacturer

Catalogue number

Authorized Representative in the European Community

Single-use only

Batch code

WARNING: THESE REAGENTS MUST NOT BE SUBSTITUTED FOR THE MANDATORY POSITIVE AND NEGATIVE SAMPLE REAGENTS PROVIDED WITH MANUFACTURED TEST KITS.

REDx™ STI Negative Control

FOR IVD USE.

INTENDED USE

REDx™ STI Negative Control is an unassayed control intended to monitor laboratory testing performance, procedures, and workflow with nucleic acid tests (NATs) that detect sexually transmitted infections (STI) (Human papillomavirus (HPV), *Mycoplasma genitalium* (MG), *Trichomonas vaginalis* (TV) and *Neisseria gonorrhoeae* (NG)) in human cervical, endo-cervical, vaginal, penile and urine samples collected in a methanol based transport medium.

PRODUCT DESCRIPTION

REDx™ STI Negative Control is formulated with human fibroblast cells. REDx™ STI Negative Control can be utilized as an external sample to monitor the process of STI nucleic acid detection assays including extraction and purification, amplification, and detection.¹

REDx™ STI Negative Control does not have an assigned value ("unassayed"). It is required that each laboratory establish an acceptance range for each lot of REDx™ STI Negative Control and for each assay procedure with which it is intended to be run, prior to routine use in the laboratory⁴.

PRINCIPLES OF THE PROCEDURE

REDx™ STI Negative Control is designed as an external independent sample for use with laboratory testing of HPV, MG, TV and NG sexually transmitted infections and provides sample adequacy nucleic acid targets, according to ISO15189 and CLIA regulations. REDx™ STI Negative Control is manufactured from human fibroblast cells. The fibroblast cells are suspended in an alcohol-based transport medium.

REAGENTS

Cat. No RED-99-M1 1 vial containing inactivated human fibroblast cells.

NOTE: REDx™ STI Negative Sample contains a suspension of cells in an alcohol solution and may therefore exhibit slight cloudiness.

LIMITATIONS OF THE PROCEDURE

REDx™ SAMPLES MUST NOT BE SUBSTITUTED FOR THE POSITIVE AND NEGATIVE SAMPLE REAGENTS PROVIDED WITH MANUFACTURED TEST KITS.

TEST PROCEDURES and INTERPRETATION OF RESULTS provided by manufacturers of test kits must be followed closely.

Deviations from procedures recommended by test kit manufacturers may produce unreliable results.

REDx™ STI Negative Control DOES NOT HAVE AN ASSIGNED VALUE and may not be suitable for use with all STI test kits and procedures. Specific levels of reactivity, when Sample Adequacy Control (SAC) is used, will vary among different manufacturers' assays, different procedures and different laboratories. Procedures for implementing a quality assurance program and monitoring test performance on a routine basis must be established by each individual laboratory. Each laboratory is required to establish its own range of acceptable values⁴.

Samples are not calibrators and should not be used for assay calibration.

REDx™ STI Negative Control is recommended for use with nucleic acid test only.

Adverse shipping and storage conditions or use of outdated samples may produce erroneous results.

REDx™ STI Negative Control might not be suitable for nucleic acid tests without an extraction step.

WARNINGS AND PRECAUTIONS

For IVD use.

For Professional and Trained Laboratory Personnel Use Only

Safety Precautions

1. Raw material used for REDx™ STI Negative Control preparation is inactivated.
2. Use Centers for Disease Control and Prevention (CDC) recommended universal precautions for handling the samples and human specimens².
3. Do not pipette by mouth; do not smoke, eat or drink in areas where specimens are being handled. Clean any spillage by immediately wiping up with 0.5% sodium hypochlorite solution. Dispose of all specimens, samples and materials used in testing as though they contain infectious agents.
4. REDx™ STI Negative Control must be disposed of by following RCRA ID#D001 guidelines for ignitable waste³.

5. Keep REDx™ STI Negative Sample closed when not in use; avoid direct inhalation of the solution and use with ventilation.

Handling Precautions

1. Do not use samples beyond the expiration date.
2. Avoid contamination of samples when opening and closing the vials
3. REDx™ STI Negative Control diluent contains a FLAMMABLE liquid; keep away from all sources of ignition.

STORAGE INSTRUCTIONS

Store REDx™ STI Negative Control at 2-8°C until use.

Once opened REDx™ STI Negative Control should not be reused. Store vials upright to prevent leakage.

MATERIALS PROVIDED

REDx™ STI Negative Control – 1 vial x 1 mL

MATERIALS REQUIRED, BUT NOT PROVIDED

Refer to the instructions supplied by manufacturer of the test kit to be used.

PROCEDURE

When including the REDx™ STI Negative Control in a test run, the exact same procedure for unknown specimens collected in STM (sample transport medium) must be used. Refer to the manufacturer supplied instructions for use provided with the STI test kit.

1. Mix by vortexing the REDx™ STI Positive Control immediately prior use.
2. Skip the swab elution step (described in some assay instructions for use)
3. Use 200-1000 µL from the sample for the kit nucleic acid extraction step.
4. After extraction proceed with RT-PCR by using the sample test volume specified in the assay instructions for use (usually 5-20 µL from the eluted purified nucleic acid volume).

NOTE: Samples must NOT be substituted for the internal kit positive and negative controls. REDx™ STI Negative Control is an independent control intended to monitor laboratory testing performance. However, in case of any serious incident that has occurred in relation to the product, report the event to Microbix and the competent authority of the country in which the event occurred.

Levels of reactivity of REDx™ STI Negative Control, when Sample Adequacy Control (SAC) is used, may vary with different manufacturer's tests and different test kit lots. As REDx™ STI Negative Control does not have an assigned value, the laboratory is required to establish an acceptance range for each lot of REDx™ STI Negative Control.

TROUBLESHOOTING

When results for the REDx™ STI Negative Control are outside of the established laboratory acceptance range for SAC, it may be an indication of unsatisfactory test performance.

Possible sources of error include: deterioration of test kit reagents, operator error, faulty performance of equipment, or contamination of reagents; internal laboratory procedures should be followed.

REFERENCES

1. *Accurate Results in the Clinical Laboratory 2013*, ISBN: 978-0-12-415783-5
2. *CDC Recommendations for prevention of HIV transmission in health care settings. MMWR 36 (suppl. 2), 1987.*
3. *Treatment standards for hazardous waste; 40 CFR 268.40 Subpart D. D001: Ignitable characteristics of waste.*
4. *Kinns H, Pitkin S, Housley D, et al. J Clin Pathol 2013;66:1027-1032.*
5. *Statistical Quality Sample for Quantitative Measurements: Principles and Definitions; Approved Guideline- Second Edition. NCCLS document C24-A2, 1999.*

For assistance, contact Microbix Technical Support at +1-905-361-8910.

P/N RED-99-M1.5R1

08 October 2024

PHONE: 905.361.8910

FAX: 905.361.8911

Toll Free (US only): 800.794.6694

www.microbix.com

CERTIFICATE OF ANALYSIS
REDx™ STI Negative Control

Product Information

| | |
|---------------------|---|
| Product Name | REDx™ STI Negative Control |
| Catalogue Number | RED-99-M1 |
| Strain/Cell line | MRC-5 |
| Lot number | 99M10018A1 |
| Analyte Value | Unassayed. Value not assigned by the manufacturer |
| Product Description | Unassayed qualitative negative control for the detection of <i>Mycoplasma genitalium</i> (MG), Human Papillomavirus (HPV), <i>Neisseria gonorrhoeae</i> (NG) and <i>Trichomonas vaginalis</i> (TV). Contains human DNA. |
| Fill Volume | 1 mL |
| Buffer | Aqueous solution containing methanol as preservative |
| Storage conditions | 2-8°C |
| Expiry Date | 01 May 2028 |
| Intended Use | For In Vitro Diagnostic Use. |

Quality Control Information

| | |
|-------------------------|---------------------|
| Test method for MG | ResistancePlus MG |
| Test method for TV | Cepheid Xpert TV |
| Test method for HPV/HBB | Cepheid Xpert HPV |
| Test method for NG | Cepheid Xpert CT/NG |
| Test Results | Observed Result |
| MG | Negative |
| TV | Negative |
| HPV | Negative |
| NG | Negative |
| hDNA (HBB) | Positive |
| Inactivation | Inactivated |

This product is unassayed control. The reported QC information is not intended to represent product specification in any commercial or lab developed test. Laboratory should establish its own analyte value and ranges.

Microbix Biosystems Inc.
265 Watline Avenue
Mississauga, Ontario, L4Z 1P3
Obelis S.A.
Bd. General Wahis, 53 1030 Brussels
Belgium

Information

Please ensure that your institution is consulted before recycling an item.

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REFERENZEN

1. *Accurate Results in the Clinical Laboratory 2013*, ISBN: 978-0-12-415783-5
2. *CDC Recommendations for prevention of HIV transmission in health care settings*. MMWR 36 (suppl. 2), 1987.
3. *Treatment standards for hazardous waste; 40 CFR 268.40 Subpart D. D001: Ignitable characteristics of waste*.
4. Kinns H, Pitkin S, Housley D, et al. *J Clin Pathol* 2013;66:1027–1032.
5. *Statistical Quality Sample for Quantitative Measurements: Principles and Definitions; Approved Guideline– Second Edition*. NCCLS document C24-A2, 1999.

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REFERENCIAS

1. *Accurate Results in the Clinical Laboratory 2013*, ISBN: 978-0-12-415783-5
2. *CDC Recommendations for prevention of HIV transmission in health care settings*. MMWR 36 (supl. 2), 1987.
3. *Treatment standards for hazardous waste; 40 CFR 268.40 Subpart D. D001: Ignitable characteristics of waste*.
4. Kinns H, Pitkin S, Housley D, et al. *J Clin Pathol* 2013;66:1027-1032.
5. *Statistical Quality Sample for Quantitative Measurements: Principles and Definitions; Approved Guideline- Second Edition*. NCCLS document C24-A2, 1999.

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BRONNEN

1. *Accurate Results in the Clinical Laboratory 2013*, ISBN: 978-0-12-415783-5
2. *CDC Recommendations for prevention of HIV transmission in health care settings*. MMWR 36 (suppl. 2), 1987.
3. *Treatment standards for hazardous waste; 40 CFR 268.40 Subpart D. D001: Ignitable characteristics of waste.*
4. Kinns H, Pitkin S, Housley D, et al. *J Clin Pathol* 2013;66:1027-1032.
5. *Statistical Quality Control for Quantitative Measurements: Principles and Definitions; Approved Guideline- Second Edition*. NCCLS document C24-A2, 1999.

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R-Biopharm AG